

# TRUPOTAN<sup>®</sup> G

## replacement syntan

Basis:	phenolic compounds
Appearance:	yellow brown powder
Charge:	anionic
pH value(1 : 10):	approx. 3
Light fastness:	good
Heat yellowing:	good

### **Properties:**

**TRUPOTAN G** is a replacement syntan for the retanning of shoe upper leather and suede.

**TRUPOTAN G** brings about an increase in the fullness, especially in the more opened structured belly areas.

**TRUPOTAN G** improves the grain tightness.

**TRUPOTAN G** levels the dyeings without bleaching.

### **Application:**

In the retanning normally 3 to 8 % **TRUPOTAN G** is recommended, based on the shaved weight. **TRUPOTAN G** encourages further penetration of the vegetable retannage, syntans and polymers.

### **Safety and storage:**

When handling **TRUPOTAN G**, normal safety precautions associated with the handling of chemicals should be observed. For more specific details please refer to our safety data sheets.

**TRUPOTAN G** can be stored for up to 24 months, if kept in a dry place.